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WORKING INSTRUCTIONS MODELCAST S

Modelcast S is extremely fine grained, temperature resistant phosphate investment material used for making partial dentures with metal base for speed or conventional heating method. It can be duplicated only in silicone. For good mixture, good chemical binding and smooth cast the investment should be mixed 15 seconds by hand and 60 seconds by vacuum.

- Before you mix the powder and liquid the mixing bowl has to be clean and dry.
- When you use the silicon for duplicating and you place it under the pressure, you must continue the complete process under pressure.
- Use different mixing bowl for investment and stone. Take care that the investment does not come into the contact with stone.
- We recommend using degreasing agent for wax parts. (Interwaxit). Interwaxit is necessary well-dry after use.
- Keep furnace door closed for 10 minutes after insertion of the mould!

WORKING INSTRUCTIONS:

- First pour liquid and then add powder.
- Manual mixing time 15 s
- Vacuum mixing time 60 s
- Duplicate model in silicone.
- Prepare the duplicated model (The base of the model should have a minimum thickness of 1 cm) –
- Remove the model from silicon after 30- 40 minutes.
- Dip/Spray the duplicated model (100 °C): MODELHART SPRAY
- Do the wax modelation
- Invest into the mould (see the mixing ratio).
- Scratch the surface of the mould before place it in the furnace

MIXING RATIO POWDER/ LIQUID (DUPLICATION IN SILICONE)

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SPEED HEATING METHOD

Place the mould in the furnace at 900- 950 °C 17 -20 minutes after first contact of liquid and powder.

Holding time minimum .60 minutes – (depending on size and number of moulds).

CONVENTIONAL HEATING METHOD:		
Temperature	Rate of heating	Holding time
Room temperature – 270 °C	5 °C/мин	40 min
270 -580 °C	7 °C/мин	30 min
580 -950 °C	9 °C/мин	60 min

Holding time minimum 60 minutes – (depending on size and number of moulds).

Casting (I-MG, I-MG FH, I-MG EKO, I-MG HE)

Let the mould cool down to the room temperature, do not force the cooling.

After cooling to the room temperature take the cast out of the investment.

STORAGE: at room temperature

INVESTING TEMPERATURE: 19 °C -23 °C

SAFETY INSTRUCTIONS: Product contains quartz. Risk of damage to lungs through prolonged or repeated exposure if inhaled. Do not breathe dust. We recommend wearing half facemask with filter FFP1.

WARRANTY! Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to further development. Therefore changes in construction and composition are reserved.

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